Product Quality

 **LABinside**

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PRODUCTS
Quality is the most important criterion in the manufacturing
of our products.

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Quality is the most important criterion in the manufacturing of our products. For this purpose, we have implemented a control system which monitors the physical, chemical and biological quality of our products. The checking of quality starts with the incoming raw material, is continued in production processes and ends with the dispatch of the finished product. Strict controls, conducted according to the legal provisions and specific standards, accompany the product. All products can be traced back to their production period, machine operators and tools, all the way back to the starting material used. With each delivery of high quality goods, our customers receive a product specific certificate, which confirm our high quality standards.

For this reason LABinside guarantees for a multitude of its products a routine control of raw materials (→p. 7). The products are checked for the following quality relevant characteristics: Free from endotoxins (→p. 13), free from DNA/DNase/RNase (→p. 9). For all our sterile products, we guarantee a SAL 10^{-3} (→p. 11).

Particularly in the field of in-vitro diagnostics toxic effects on cells or tissue cultures plays an important role. It is vital therefore that any toxic effect during the cultivation and storage of cells or tissues for in-vivo use are eliminated(→p. 15).

This service will be performed to ensure your security and is free of charge. Each product tested will be delivered with a certificate of quality, like that shown on page 5.

→To order the detailed up-to-date listing of products tested for cytotoxicity, please give us a call or send an email.

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Control of raw materials	p. 6 - 7
Free from DNA	p. 8 - 9
Free from DNase	p. 8 - 9
Free from RNase	p. 8 - 9
Sterility SAL 10^{-3}	p. 10 - 11
Free from endotoxins	p. 12 - 13
Free from cytotoxicity	p. 14 - 15

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QUALITY
Quality assurance marks not only guarantee the quality but also give
a guideline for the products.

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PLASTIC GRANULES
The granules are homogenised at 240 to 290°C and
then further processed to high quality products by
injection moulding.

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
LABinside products are manufactured from raw materials that:

- meet the requirements of the FDA directive 21CFR 177.1640
“Polystyrene and rubber-modified polystyrene”
- meet the EU directive 94/62/EC, this regulation does require that the maximum level for heavy metals is below 100 ppm.

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HUMAN HAIRS
Contamination with biological material will generate
false results in DNA-analysis.

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DNase/RNase

DNases (deoxynucleases) and RNases (ribonucleases) are enzymes that hydrolyse DNA (deoxyribonucleic acid) and RNA (ribonucleic acid) into small fragments, thus destroying them. The presence of these enzymes in laboratory products used for PCR (polymerase chain reaction) can lead to false results due to the reduction or absence of in-vitro transcription reactions.

Detection of these enzymes is conducted by PCR:

The products are washed with 0.5 % Tween20 to completely detach enzymes that might be adhering to any surfaces. The wash solution is incubated together with a DNA (or RNA) standard and the residual DNA (or RNA) content is determined by PCR, or RT-PCR (for RNase) detection. A positive control (DNA or RNA standard without wash solution) and a negative control (H2O, 0.5 % Tween20 and buffer used) are included in each test. The amplification products are detected online by means of FRET (fluorescence resonance energy transfer).

(Human) DNA

Deoxyribonucleic acid is a biopolymer and the carrier of genetic information. For accurate, meaningful experimental results when using techniques such as PCR, DNA sequencing or other molecular biology techniques, it is vital to avoid contamination with foreign DNA.

Detection of human DNA is conducted by PCR:

The products are washed with 0.5 % Tween20 to completely detach nucleotides that may be adhering to any surfaces. The wash solution is tested for the presence of DNA by means of PCR. A positive control (DNA standard) and a negative control (0.5 % Tween20) are included in each test. The amplification products are detected online by means of FRET.

Ribonucleases and DNA are found in all biological materials (skin, hair, etc.) and both the enzymes and the DNA are extremely resistant to environmental influences. Sterilisation of products by irradiation will only partially inactivate or destroy them. Experimental evidence shows that even doses of 20 kGy and above only induce low levels of inactivation in these molecules. Therefore, it is vital to prevent contamination due to human or animal material.

Contamination with nucleases and DNA due to contact with human and animal material must therefore be prevented.

Our production environment is designed in such a way as to minimise the risk of contamination.

- All production staff wear overalls, gloves and a hairnet, in order to prevent direct contact of the products with human material.
- The ambient air in the production area is strictly controlled by appropriate measures and is monitored for particles and microorganisms.

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MICROORGANISM
Sterile conditions are essential for cell and tissue cultivation.

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Sterile means free of viable microorganisms.
 Sterilisation is performed to inactivate viable microorganisms.
 However, the inactivation of microorganisms by sterilisation often
 approximates to an exponential function. Therefore, the presence of viable
 microorganisms is defined as a “probability”. This probability can be
 reduced to a minimum by sterilisation, but it can never be equal to zero.
 In the sterilisation procedure, the probability of a non-sterile item is
 referred to as the “sterility assurance level” SAL.

For our sterile products, we guarantee a SAL10⁻⁶.
 The radiation dose is determined in accordance with ISO 11737.
 For calculation of the radiation dose required to guarantee this SAL10⁻⁶.
 it is necessary to determine the population of viable microorganisms
 (bioburden) prior to sterilisation. For this purpose, a defined number of
 items are taken from three production units to determine the bioburden.
 However, since some microorganisms can adhere to the surface of the
 items a so-called recovery rate have to be evaluated advance for each
 individual product. For this purpose, sterile products are contaminated
 with a defined number of microorganisms and then the bioburden rate is
 determined. The correction factor thus determined and the bioburden of
 the non-sterile products are used to calculate the radiation dose.
 After validation of the radiation dose determined, the products are
 irradiated accordingly. The bioburden for each product is checked at
 regular intervals and re-evaluated whenever any relevant changes are
 made to the production process.

Sterilisation is monitored with dosimeters and validated at regular intervals
 in accordance with ISO 11137, EN 552, ANSI/AAMI ST31/ST32.

Each sterilisation is accompanied by documentation of the radiation dose
 measured.

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For our sterile products we guarantee a SAL 10⁻⁶.

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HORSESHOE CRAB
Endotoxins are detected unequivocally by means of the blue blood of the horseshoe crab (*limulus polyphemus*).

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Endotoxins are a subset of pyrogens that originate from gram-negative bacteria. They are a lipopolysaccharide complex occurring in the outer layer of the bilayered gram-negative bacterial cell, and can include a wide variety of different biological reactions in humans, ranging from fever to septic shock.

Endotoxins are extremely resistant to environmental influences such as heat and are not affected by simple heat-sterilisation or drying. They can only be destroyed by dry heat at over 200 ° C, with a incubation time of at least one hour.

The detection of endotoxins is conducted by means of the Limulus amoebocyte lysate test (LAL). The test utilises the blood cells from the horseshoe crab (Limulus polyphemus) which react with endotoxins to produce a gelling of the lysate. This phenomenon confirms the presence of endotoxins in the sample and standard curves are used to determine the sensitivity of the test system.

All products tested are free from endotoxins up to a detection limit of ≤ 0.06 EU/ml (EU=endotoxin units).

The test system meets the FDA guidelines for medical products (12/87).

The U.S. Pharmacopoeia stipulates that all medical products must not exceed the limit value of 0.5 EU/ml. Products that come into contact with cerebrospinal fluid must not exceed ≤ 0.06 EU/ml.


Contamination with endotoxins due to contact with human and animal material must therefore be prevented.

Our production environment is designed in such a way as to minimise the risk of contamination:

- All production staff wear overalls, gloves and a hairnet, in order to prevent direct contact of the products with human material.
- The ambient air in the production area is strictly controlled by appropriate measures and is monitored for particles and microorganisms.

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All products tested are free from endotoxins up to a detection limit of ≤ 0.06 EU/ml.

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HEPATOCTYES
A cytotoxic environment results in cell necrosis.

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BHK-21 cell은 MEM(Earle's salts, 10% FBS, 4 mM Glutamine, 1 x MEM non-essential amino acids)배지로 ,
Balb/3T3 cell은 Dulbecco's MEM(10% FBS, 4 mM Glutamine)배지로,

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- Test samples: 1.25 cm² sample material/ml extraction medium
Controls: Negative control: 1.25 cm² control material/ml extraction
 medium, HDPE from U.S. Pharmacopoeia
 Positive control: 1.25cm² control material/ml extraction
 medium, Organo-tin PVC from Portex
 Alternative: 0.1% phenol in the extraction medium
 Reagent control: cell culture medium
 Extraction time: 72h at 37°C and 5% CO₂.

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In many experiments toxic effects on cells or tissue cultures play an important role particularly in the field of in-vitro diagnostics. It is vital therefore that any toxic effects during the cultivation and storage of cells or tissues for in-vivo use are eliminated.

The relevant culture vessels and all objects that come into direct contact with the cells and tissues must therefore be free of cytotoxic substances.

The detection of cytotoxicity is evaluated with mammalian cells according to EN ISO 10993.

Any adverse biological reaction due to contact with our production materials e.g. polystyrene or polypropylene is closely monitored by culturing cells in the presence of extracts from these polymers.

The following cells are used:

BHK-21 cells (syrian hamster kidney; DSMZ ACC 61) are cultivated as follows:

Medium . MEM with Earle's salts with 10% FBS, 4mM Glutamine and 1 x MEM non-essential amino acids.

Balb/3T3 cells (mouse, embryo; ATCC CCL163) are cultivated as follows:

Medium . Dulbecco's MEM with 10% FBS, 4mM Glutamine.

The cells are cultivated in non-cytotoxic cell culture flasks at 37°C in 5% CO₂.

The materials to be tested are prepared as extracts under sterile conditions, according to ISO 10993-12:

- Test samples: 1.25 cm² sample material/ml extraction medium
Controls: Negative control: 1.25 cm² control material/ml extraction medium,
 HDPE from U.S. Pharmacopoeia
 Positive control: 1.25cm² control material/ml extraction medium,
 Organo-tin PVC from Portex
 Alternative: 0.1% phenol in the extraction medium
 Reagent control: cell culture medium
 Extraction time: 72h at 37°C and 5% CO₂.

Test for cytotoxicity of the extracts (according to ISO 10993-5):

The cell culture medium of the cultivated cells is exchanged with the same amount of the extracts and incubated for at least 24 hours at 37°C and 5% CO₂. After the incubation period, all of the cells are checked for their morphology and cell count:

Laboratory T.O.T.A.L Solution

www.LABinside.com